

EASY PULSE

Clinical Evidence – Study Background Information

“Head-up position significantly improves CPR”

According to a recent study a tilt angle during cardiopulmonary resuscitation improves brain flow and is more effective.^{1,2}

During CPR, there is diminished venous drainage due to the body being level. Venous drainage occurs passively, so it is improved by raising the head above the heart. CPR with head-up positioning provides a new way to augment cerebral perfusion without impairing perfusion of the heart.

SCHILLER EASY PULSE meets this challenge: **SCHILLER EASY PULSE is the only device that can perform CPR in any position.** The device is strapped to the body and on a backboard, which allows following the body's movement even in extreme conditions, for example during air rescue or with a patient trapped in a car seat after an accident.



1 Resuscitation 87 (2015) 38-43

Tilting for perfusion: Head-up position during cardiopulmonary resuscitation improves brain flow in a porcine model of cardiac arrest

Guillaume Debaty^{a,b}, Sang Do Shin^c, Anja Metzger^e, Taeyun Kim^c, Hyun Ho Ryu^d, Jennifer Rees^a, Scott McKnite^e, Timothy Matsuura^f, Michael Lick^e, Demetris Yannopoulos^a, Keith Lurie^e

^a Department of Medicine-Cardiovascular Division, University of Minnesota, Minneapolis MN, USA

^b UJF-Grenoble 1/CNRS/CHU de Grenoble/TIMC-IMAG, UMR 5525, Grenoble, France

^c Department of Emergency Medicine, Seoul National Hospital, Seoul, Republic of Korea

^d Department of Emergency Medicine, Chonnam National University Hospital, Republic of Korea

^e Department of Emergency Medicine, Hennepin County Medical Center, University of Minnesota, Minneapolis, MN, USA

^f Department of Integrated Biology & Physiology, Medical School, University of Minnesota, Minneapolis, MN, USA

2 JEMS Journal of Emergency Medical Services

Tilt Angle Significantly Affects CPR

Sat, Mar 7, 2015 1 By Keith Wesley, MD, FACEP, Karen Wesley, NREMT-P []

Debaty G, Shin SD, Metzger A, et al. Tilting for perfusion: Head-up position during cardiopulmonary resuscitation improves brain flow in a porcine model of cardiac arrest. Resuscitation. 2015;87:38-34.



Americas
SCHILLER Americas Inc.
 Doral, Florida 33172

North America:
 Phone +1 786 845 06 20
 Fax +1 786 845 06 02
sales@schilleramerica.com
www.schilleramerica.com

Latin America & Caribbean
 Phone +1 305 591 11 21
 Fax +1 786 845 06 02
sales@schillerla.ch
www.schillerla.ch



Asia
SCHILLER Asia-Pacific / Malaysia
 52200 Kuala Lumpur, Malaysia
 Phone +603 6272 3033
 Fax +603 6272 2030
sales@schiller.com.my
www.schiller-asia.com



Austria
SCHILLER Handelsgesellschaft m.b.H.
 A-4040 Linz
 Phone +43 732 70 99
 Fax +43 732 757 000
sales@schiller.at
www.schiller.at



China
Alfred Schiller (Beijing) Medical Equipment Co., Ltd.
 100015 Beijing, China
 Phone +86 010 52007020
 Fax +86 010 52007020-8016
info@schillerchina.com
www.schillermedical.cn



Croatia
SCHILLER d.o.o.
 10000 Zagreb
 Phone +385 1 309 66 59
 Fax +385 1 309 66 60
info@schillerzg.hr
www.schiller.hr



France
SCHILLER Médical S.A.S.
 F-67162 Wissembourg
 Phone +33 3 88 63 36 00
 Fax +33 3 88 63 36 49
info@schiller.fr
www.schiller-medical.com



France (distribution France)
SCHILLER France S.A.S.
 F-77608 Bussy St Georges
 Phone +33 1 64 66 50 00
 Fax +33 1 64 66 50 10
contact@schillerfrance.fr
www.schiller-france.com



Germany
SCHILLER Medizintechnik GmbH
 D-85622 Feldkirchen b. München
 Phone +49 89 62 99 81 0
 Fax +49 89 609 50 90
info@schillermed.de
www.schillermed.de



Hungary
SCHILLER Diamed Ltd.
 H-1141 Budapest
 Phone +36 1 383 4780 / +36 1 460 9491
 Fax +36 1 383 4778
info@schillerhungary.hu
www.schillerhungary.hu



India
SCHILLER Healthcare India Pvt. Ltd.
 Mumbai - 400 059, India
 Phone +91 22 6152 3333 / 2920 9141
 Fax +91 22 2920 9142
sales@schillerindia.com
www.schillerindia.com



Poland
SCHILLER Poland Sp. z o.o.
 PL-02-729 Warszawa
 Phone +48 22 843 20 89 / +48 22 647 35 90
 Fax +48 22 843 20 89
schiller@schiller.pl
www.schiller.pl



Russia & C.I.S.
AO SCHILLER.RU
 125124 Moscow, Russia
 Phone +7 (495) 970 11 33
 Fax +7 (495) 956 29 10
mail@schiller.ru
www.schiller.ru
www.schiller-cis.com



Serbia
SCHILLER d.o.o.
 11010 Beograd
 Phone +381 11 39 79 508
 Fax +381 11 39 79 518
info@schiller.rs
www.schiller.rs



Slovenia
SCHILLER d.o.o.
 2310 Slovenska Bistrica
 Phone +386 2 843 00 56
 Fax +386 2 843 00 57
info@schiller.si
www.schiller.si



Spain
SCHILLER ESPAÑA, S.A.
 ES-28232 Las Rozas/Madrid
 Phone +34 91 713 01 76
 Fax +34 91 355 79 33
schiller@schiller.es
www.schiller.es



Switzerland
SCHILLER-Reomed AG
 CH-8953 Dietikon
 Phone +41 44 744 30 00
 Fax +41 44 740 37 10
sales@schiller-reomed.ch
www.schiller-reomed.ch



Turkey
SCHILLER TÜRKİYE
 Okmeydanı-Sisli - Istanbul
 Phone +90 212 210 8681 (pbx)
 Fax +90 212 210 8684
info@schiller.com.tr
www.schiller-turkiye.com



United Kingdom
SCHILLER UK
 Bells Hill, ML4 3PR
 Phone +44 1698 744 505
 Fax +44 1698 744 474
sales@schilleruk.com
www.schilleruk.com



Suomi - Finland
Medifin Healthcare Oy
 Puh. +358 (0)10 321 3870
asiakaspalvelu@medifin.fi
www.medifin.fi



SCHILLER
 The Art of Diagnostics



Tilt Angle Significantly Affects CPR

Sat, Mar 7, 2015 | By [Keith Wesley, MD, FACEP](#) , [Karen Wesley, NREMT-P](#) []

The Research

Debaty G, Shin SD, Metzger A, et al. Tilting for perfusion: Head-up position during cardiopulmonary resuscitation improves brain flow in a porcine model of cardiac arrest. *Resuscitation*. 2015;87:38-34.

The Science

The authors of this study had previously noted elevated pressures inside the brain, intracranial pressure (ICP), and hypothesized that raising the head during resuscitation could lower said pressure. They studied 22 pigs that were placed into v fib. The pigs and a LUCAS mechanical CPR device were bolted to a backboard and secured to a table that could be tilted in increments from 0-50 degrees. All pigs were placed in untreated v fib for six minutes followed by two minutes of CPR with an impedance threshold device (ITD) in a horizontal position.

After six minutes of untreated v fib, CPR was performed on 14 pigs at 0 degrees, 30 degrees head up and 30 degrees head down. Microspheres were used to measure organ blood flow in eight pigs. ITD CPR was performed on eight additional pigs at 0, 20, 30, 40 and 50 degrees head up.

They discovered that coronary and cerebral perfusion increased significantly when tilting the table upwards. With 0, 10, 20, 30, 40 and 50 degrees head up tilt, ICP values were 21 ± 2 , 16 ± 2 , 10 ± 2 , 5 ± 2 , 0 ± 2 , -5 ± 2 respectively ($p <$

0.001). Tilting the head down resulted in significantly elevated ICP and decreased brain perfusion.

They concluded that cerebral perfusion was significantly affected by the tilt angle of the pig during CPR.

Doc Wesley Comments

My first reaction when I read this study was, "Why didn't I think of that?" We've known for decades that raising the head of the bed is the first step in lowering the pressure in the brain of patients with head injuries. Cardiac arrest research was initially oriented with the goal of resuscitating the heart. Today, we realize there's no benefit to getting the heart back if the patient is left with poor neurological outcome.

Blood flow to the brain occurs during systole, or in an arrest situation, when we compress on the chest. During CPR, there's diminished venous drainage due to the body being level. Venous drainage occurs passively but is promoted by raising the head above the heart. The placement of an ITD enhances venous drainage and promotes blood return to the heart. Each time we compress on the chest, a hammerhead of pressure pounds the brain—which is already in shock. This has led some to theorize that traditional CPR may actually harm the brain, as a result of increased cerebral pressure and decreased perfusion.

In this study, the researchers removed the LUCAS for a short period during the 30-degree head up phase and documented that it significantly augmented the effect of the tilt. These same researchers have also completed a study showing that tilting the head upward improves neurological outcome in resuscitated pigs.

Earlier this year, Karen and I, along with several medical directors from across the country, had the opportunity to witness this effect first-hand at an event hosted by Advanced Circulatory System, the maker of the ResQPOD ITD (of which I have no financial involvement).

Frankly? I was blown away. It just makes sense. Other medical directors and I are in the process of determining how it could be implemented in the field. Using current technology, it seems quite possible to tilt the body at least 30 degrees during both manual and mechanical CPR. This could very well become the new standard for CPR delivery.

Medic Wesley Comments


After 34 years in EMS, and multiple changes in CPR and resuscitation guidelines, I tend to take new updates with a "ho-hum" attitude. This time, I was able to witness the science.

By tilting the head of the subject during compression, the increased ICP seen with manual and mechanical CPR wasn't present. It's so simple, as Doc states. Why hasn't it been thought of before?

We seemed so focused on timing of breaths, quality compressions and early defibrillation that the mechanical damage done to the brain never seemed to come up. Pushing hard and fast may help the heart but at the same time it may actually be harming our patient's brain. Now that's just a supposition on my part. But it does bode well for looking at everything from—pardon the pun—a different angle.

Using the current science of compressions, along with the new information suggested in this study, this may be the complete package when it comes to successful resuscitation with good neurological outcome. Any of you who are inventors should probably start figuring out a device for CPR that supports the mechanical device and elevates the patient's head. After seeing the success of this study, I'm relatively certain we'll be looking for a way to implement it in our ambulances soon.

By

**Keith Wesley, MD, FACEP**

Keith Wesley is a board-certified emergency medicine physician living in Wisconsin. Originally from Tyler, Texas, he graduated from Brigham Young University in 1982 and Baylor College of Medicine in Houston in 1986. He completed an Emergency Medicine Residency at Methodist Hospital in Indianapolis, where he gained his first exposure to EMS flying air medical missions. Dr. Wesley has been involved in EMS since 1989 working with many services in Wisconsin. In 1992 he was selected by the Governor as a founding member of the Wisconsin State Physician Advisory Committee and served for 12 years, the past four years as chair. In 2006, Dr. Wesley was selected as the Wisconsin State EMS medical director and continues to provide medical oversight to several services throughout Wisconsin. He is also the past Chair of the National Council of State EMS Medical Directors and is active in National Association of EMS Physicians. In 2007 Dr. Wesley became the Minnesota State EMS medical director and medical director for HealthEast Medical Transportation in St. Paul. Dr. Wesley has authored four EMS textbooks and numerous articles and papers on EMS and is a frequent speaker at state and national EMS conferences. Dr. Wesley is the owner of EMS Consulting and Education, which provides on-line educational programs for EMS providers and medical directors, as well as consultative services for EMS delivery and quality improvement.

**Karen Wesley, NREMT-P**

Karen Wesley, NREMT-P, is a paramedic and educator for Mayo Clinic Medical Transport and is the medic team leader for the Eau Claire County (Wis.) Regional SWAT team. She can be reached at admKaren22@hotmail.com.

Copyright © 2017: PennWell Corporation, Tulsa, OK. All Rights Reserved.



ELSEVIER

Contents lists available at ScienceDirect

Resuscitation

journal homepage: www.elsevier.com/locate/resuscitation

Experimental paper

Tilting for perfusion: Head-up position during cardiopulmonary resuscitation improves brain flow in a porcine model of cardiac arrest[☆]



Guillaume Debaty^{a,b,*}, Sang Do Shin^c, Anja Metzger^e, Taeyun Kim^c, Hyun Ho Ryu^d, Jennifer Rees^a, Scott McKnite^e, Timothy Matsuura^f, Michael Lick^e, Demetris Yannopoulos^a, Keith Lurie^e

^a Department of Medicine-Cardiovascular Division, University of Minnesota, Minneapolis, MN, USA

^b UJF-Grenoble 1/CNRS/CHU de Grenoble/TIMC-IMAG UMR 5525, Grenoble, France

^c Department of Emergency Medicine, Seoul National Hospital, Seoul, Republic of Korea

^d Department of Emergency Medicine, Chonnam National University Hospital, Republic of Korea

^e Department of Emergency Medicine, Hennepin County Medical Center, University of Minnesota, Minneapolis, MN, USA

^f Department of Integrated Biology & Physiology, Medical School, University of Minnesota, Minneapolis, MN, USA

ARTICLE INFO

Article history:

Received 22 October 2014

Received in revised form

17 November 2014

Accepted 21 November 2014

Keywords:

Cardiac arrest

Cardiopulmonary resuscitation

Mechanical CPR

Impedance threshold device

Cerebral perfusion

ABSTRACT

Introduction: Cerebral perfusion is compromised during cardiopulmonary resuscitation (CPR). We hypothesized that beneficial effects of gravity on the venous circulation during CPR performed in the head-up tilt (HUT) position would improve cerebral perfusion compared with supine or head-down tilt (HDT).

Methods: Twenty-two pigs were sedated, intubated, anesthetized, paralyzed and placed on a tilt table. After 6 min of untreated ventricular fibrillation (VF) CPR was performed on 14 pigs for 3 min with an automated CPR device called LUCAS (L) plus an impedance threshold device (ITD), followed by 5 min of L-CPR + ITD at 0° supine, 5 min at 30° HUT, and then 5 min at 30° HDT. Microspheres were used to measure organ blood flow in 8 pigs. L-CPR + ITD was performed on 8 additional pigs at 0°, 20°, 30°, 40°, and 50° HUT.

Results: Coronary perfusion pressure was 19 ± 2 mmHg at 0° vs. 30 ± 3 at 30° HUT ($p < 0.001$) and 10 ± 3 at 30° HDT ($p < 0.001$). Cerebral perfusion pressure was 19 ± 3 at 0° vs. 35 ± 3 at 30° HUT ($p < 0.001$) and 4 ± 4 at 30° HDT ($p < 0.001$). Brain–blood flow was 0.19 ± 0.04 ml min⁻¹ g⁻¹ at 0° vs. 0.27 ± 0.04 at 30° HUT ($p = 0.01$) and 0.14 ± 0.06 at 30° HDT ($p = 0.16$). Heart blood flow was not significantly different between interventions. With 0, 10, 20, 30, 40 and 50° HUT, ICP values were 21 ± 2 , 16 ± 2 , 10 ± 2 , 5 ± 2 , 0 ± 2 , -5 ± 2 respectively, ($p < 0.001$), CerPP increased linearly ($p = 0.001$), and CPP remained constant.

Conclusion: During CPR, HDT decreased brain flow whereas HUT significantly lowered ICP and improved cerebral perfusion. Further studies are warranted to explore this new resuscitation concept.

© 2014 Elsevier Ireland Ltd. All rights reserved.

1. Introduction

Generation of adequate non-invasive circulatory support during CPR remains a complex physiological challenge, with competing demands to provide sufficient venous return to refill the heart after

each compression and enough circulation to the brain to preserve neurological function.^{1–3} Blood flow is highly dependent on the refilling of the heart during the decompression phase and the vascular resistance.

Closed-chest standard (S) manual cardiopulmonary resuscitation (CPR) has traditionally been performed in the 0° supine position. Stimulated by the need to develop better techniques to transport patients in high-rise apartment buildings in small elevators with ongoing CPR in Korea, an international research collaborative was established to examine the impact of CPR in more vertical positions. Little is known about the potential benefit of tilting the head upward or downward during CPR, especially when

[☆] A Spanish translated version of the summary of this article appears as Appendix in the final online version at <http://dx.doi.org/10.1016/j.resuscitation.2014.11.019>.

* Corresponding author at: SAMU 38, Pôle Urgences Médecine Aigue, CHU de Grenoble, CS 10217, 38043 Grenoble Cedex 09, France.

E-mail address: gdebaty@gmail.com (G. Debaty).

using CPR mechanical devices known to enhance circulation and provide high quality CPR.^{4,5} We hypothesized that use of mechanical adjuncts known to enhance circulation during CPR would improve outcomes when the head and body are tilted upward vs. supine CPR or with feet up and head down.

2. Method

This study was approved by the Institutional Animal Care Committee of the Minneapolis Medical Research Foundation of Hennepin County Medical Center. All animal care was compliant with the National Research Council's 1996 Guidelines for the Care and Use of Laboratory Animals. All studies were performed by a qualified, experienced research team in Yorkshire female farm bred pigs weighing 39.3 ± 0.5 kg. A certified and licensed veterinarian assured the protocols were performed in accordance with the National Research Council's Guidelines.

2.1. Preparatory phase

The surgical preparation, anesthesia, data monitoring, and recording procedures used in this study have been previously described.^{3,6,7} Under aseptic surgical conditions, initial sedation was achieved with intramuscular ketamine (10 mL of 100 mg/mL) followed by inhaled isoflurane at a dose of 0.8–1.2%. Pigs were intubated with a 7.0 French endotracheal tube. The animal's temperature was maintained between 36.5 and 37.5°C with a warming blanket (Bair Hugger, Augustine Medical, Eden Prairie, MN). Central aortic blood pressure was recorded continuously with a electronic-tipped catheter (Mikro-Tip Transducer, Millar Instruments, Houston, TX) placed in the descending thoracic aorta. A second Millar catheter was inserted in the right atrium via the right external jugular vein. An ultrasound flow probe (Transonic 420 series multichannel, Transonic Systems, Ithaca, NY) was placed in the left common carotid artery to measure carotid blood flow (mL min^{-1}). After creating a burr hole, a Millar catheter was then inserted into the parietal lobe to measure intracranial pressure (ICP). In pigs used for the microsphere studies (see below), a second femoral artery cannulation was performed and a 7F pigtail catheter was positioned in the left ventricle under fluoroscopic guidance. All animals received an intravenous heparin bolus (100 units/kg). Animals were fasted overnight and received normal saline solution to maintain the mean right atrial pressure between 3 and 5 mmHg. The animals were ventilated with room air, using an anesthesia machine (Narkomed, Telford, PA), with a tidal volume of 10 mL/kg and a respiratory rate adjusted to continually maintain an end tidal CO_2 (ETCO_2) of 40 mmHg and O_2 saturation of >92%. Arterial blood gases (Gem 3000, Instrumentation Laboratory) were obtained at baseline, and 3 min after each change of CPR position. Surface electrocardiographic tracings were continuously recorded. All hemodynamic data including aortic pressure, right atrial pressure, ETCO_2 , ICP, and carotid blood flow were continuously monitored and recorded with a digital recording system (BIOPAC MP 150, BIOPAC Systems, Inc., CA, USA). Coronary perfusion pressure (CPP) was calculated as the difference between aortic pressure and right atrial pressure during the CPR decompression phase.⁸ Cerebral perfusion pressure (CerPP) was calculated as the difference between mean aortic pressure and mean ICP. Ultrasound derived carotid blood flow velocity was reported in mL min^{-1} . ETCO_2 , tidal volume, minute ventilation, and blood oxygen saturation were continuously measured with a respiratory monitor (COSMO Plus, Novamatrix Medical Systems, Wallingford, CT).

After the surgical preparation was complete, oxygen saturation on room air was greater than 92%, and ETCO_2 was stable between

35 and 42 mmHg for 5 min, VF was induced by delivering direct intra-cardiac current via a temporary pacing wire positioned in the right ventricle. Mechanical CPR was performed using a LUCAS 1™ (Physio-Control, Redmond, WA) compression system at a rate of 100 compressions min^{-1} with a 50% duty cycle. The LUCAS backboard was bolted to a stretcher (Stryker Corporation, Kalamazoo MI) and the pig was tied by its legs to the stretcher as well. The stretcher was attached to a tilt table built to perform CPR with different study angles. In this way the pig, stretcher, and LUCAS could be moved simultaneously while L-CPR was ongoing. An impedance threshold device with a resistance of 16 cmH_2O (ITD-16, ResQPOD™, Advanced Circulatory Systems, Roseville, MN) was attached to the endotracheal tube. Asynchronous positive pressure ventilations with supplemental oxygen at a flow of 10 L min^{-1} were delivered with a manual resuscitator bag. The tidal volume was maintained at ~ 10 mL/kg and the respiratory rate was 10 breaths min^{-1} . In addition, prior to inducing VF succinylcholine ($93.3 \text{ mcg kg}^{-1} \text{ min}^{-1}$) was administered intravenously to prevent spontaneous gasping during CPR.

Angle positions were confirmed after each change of position with a digital protractor (Mitutoyo Pro 360).

2.1.1. Protocol A

Hemodynamics and calculated coronary and cerebral perfusion pressures were the focus of Protocol A. After 6 min of untreated VF, CPR was initiated on 14 pigs with L-CPR + ITD in a 0° supine position for 3 min. This interval provided time for the hemodynamic parameters to stabilize after reperfusion. L-CPR + ITD continued thereafter without interruption for multiple sequential interventions as follows: 5 min epochs at 0°, 30° head up, and 30° head down position, an additional 2 min of L-CPR + ITD in the 30° head up position and then L-CPR alone, without the ITD, for 2 additional min while still in the 30° head up position. Pigs were then placed in the 0° supine position and defibrillated with up to three 275 joule biphasic shocks (Lifepak 15, Physio-control, Redmond, WA). Animals were then sacrificed with a 10 mL injection of saturated potassium chloride.

2.1.2. Protocol B

Cerebral and myocardial perfusion were assessed under different experimental CPR positions in Protocol B. Blood flow to the heart and brain was measured with microspheres injected into the left ventricle under baseline pre-VF conditions and during CPR as previously described.^{9–11} In the current study 15 μm diameter neutron activated Lanthanum (^{140}La), Gold (^{198}Au), Ytterbium (^{175}Yb), and Lutetium (^{177}Lu) microspheres (STERIspheres™, BioPAL™; BioPhysics Assay Laboratory, Worcester, MA) were used. Microspheres were randomly assigned for each of the respective four interventions with one type of microsphere per intervention. Microspheres were injected into a total of 8 pigs during CPR under different experimental CPR positions. The microspheres were first injected into the left ventricle under stable baseline conditions 5 min prior to the induction of VF. Then, following 6 min of untreated VF CPR was performed continuously with L-CPR + ITD for the different time intervals and using the 3 different CPR positions as described in Fig. 1. After 4 min of CPR a second microsphere was injected while the pig remained in the 0° supine position. The pig was then tilted upwards to the 30° head up position and 1 min later a third microsphere was injected. After 4 min the pig was tilted downwards in the 30° head down position and 1 min later the last microsphere was injected. The number of microspheres injected for each intervention was computed as follow:

$$\mu = 1.2 \cdot 10^6 + ((1.9 \cdot 10^5) \cdot \omega)$$




Fig. 1. Study protocol. VF: ventricular fibrillation, HUT: Head up tilt, HDT: Head down tilt

where μ = required number of microspheres; ω = pig weight.

$$\text{Baseline injection volume} = \frac{\mu}{(5 \cdot 10^8)/20}$$

$$\text{Intervention injection volume} = \frac{\mu}{((5 \cdot 10^8)/20) \cdot (5/3)}$$

Concurrently with the microsphere injections, reference blood samples were withdrawn continuously from the descending aorta at a collection rate of 10 ml min⁻¹. At the end of the procedure animals were sacrificed with potassium chloride as described above and then tissue samples from the brain (pons portion of the brain-stem, hippocampus, left and right cortex) and heart (left ventricular apex, papillary muscle, free wall) samples were obtained. Samples were desiccated and sent to the reference BioPhysics Assay Laboratory for analysis.¹¹ Data for each organ, using pooled data from the portions sampled, were reported in Section 3.

2.1.3. Protocol C

This protocol was used to determine the effect of incremental increases in the head-up angle on key hemodynamic measures. The angle-response relationship for the different hemodynamic parameters assessed was determined in 8 additional pigs as follows: after 6 min of untreated VF, L-CPR+ITD was initiated and following a 3 min stabilization period in the 0° supine position hemodynamic measurements were recorded for 1 min intervals at 0°, 10°, 20°, 30°, 40°, 50° head up tilt position. Animals were then sacrificed as described above.

2.2. Statistical analysis

Data are expressed as mean ± standard error of mean (SEM). For the primary hypotheses, mean values were compared using a Student's paired *t* test with the 0° supine reference position.

One-Way Repeated Measures Analysis of Variance with linear trend was used to compare continuous data within the different angles in Protocol C. All statistical tests were two-sided, and a *p* value of less than 0.05 was required to reject the null hypothesis. Statistical analysis was performed using IBM SPSS Statistics 21.

3. Results

3.1. Protocol A: hemodynamic parameters and perfusion pressures

Results from Protocol A are shown in Table 1 and supplementary Fig. 1. During CPR in the 30° head up position mean ICP values were nearly immediately reduced by 75% compared with the 0° supine position (7 ± 1 vs. 28 ± 2, *p* < 0.001). By contrast, the mean aortic pressure decreased by about 17% in the 30° head up position (42 ± 4 vs. 48 ± 4, respectively, *p* < 0.001). Consequently, the CerPP was significantly higher with 30° head up tilt (35 ± 3 vs. 19 ± 3, respectively, *p* < 0.001). Placement of the pig in the 30° head up position resulted in 140% reduction in decompression phase right atrial pressure relative to the 0° supine position (−4.3 ± 1 vs. 10.1 ± 1, respectively, *p* < 0.001) and an increase in the calculated CPP, as shown in Table 1 and supplementary Fig. 1. By contrast, with the 30° head down position ICP was significantly higher, while CPP and CerPP were significantly lower compared with the 0° supine position (42 ± 2 vs. 28 ± 2, 10 ± 3 vs. 19 ± 2, 4 ± 4 vs. 19 ± 3, respectively, *p* < 0.001 for each comparison). Representative pressure curves recorded during the three different positions are shown in Fig. 2. A total of 9/14 pigs were successfully defibrillated at the end of the protocol.

3.2. Protocol B: cerebral and cardiac blood flow

Protocol B was used to assess instantaneous blood flow during three L-CPR+ITD positions. Cerebral blood flow (ml min⁻¹ g⁻¹ of

Table 1
Hemodynamic parameters at baseline and between 0, 30° head up and 30° head down L-CPR+ITD.

	Baseline	CPR angle		
		0°	30°	−30°
SBP	106 ± 4	94 ± 6	83 ± 5*	90 ± 6
DBP	71 ± 4*	30 ± 3	26 ± 3*	27 ± 4
RA max	8 ± 1*	157 ± 16	123 ± 14*	172 ± 13*
RA min	0 ± 1*	10 ± 1	−4 ± 1*	17 ± 1*
CBF max	613 ± 32*	377 ± 41	329 ± 43*	327 ± 46*
CBF min	223 ± 31*	−138 ± 21	−132 ± 32	−131 ± 32
ICP max	23 ± 1*	45 ± 5	15 ± 2*	65 ± 4*
ICP min	19 ± 1*	15 ± 2	−2 ± 2*	24 ± 1*

SBP: systolic blood pressure, DBP: diastolic blood pressure, RA: right atrial pressure, CBF: carotid blood flow; ICP: intracranial pressure (maximum and minimum value during compression and decompression).

* *p* < 0.05 compared with 0° supine CPR.




Fig. 2. Representative pressure curve of aortic pressure (AO), right atrial pressure (RA), intracranial pressure (ICP), coronary perfusion pressure (CPP) and cerebral perfusion pressure (CerPP) during 0° CPR, 30° Head up CPR and 30° head down CPR (−30°).

tissue) during 0° supine L-CPR + ITD was 35% of baseline (0.19 ± 0.04 vs. 0.54 ± 0.07 , $p = 0.03$) and cardiac blood flow was 28% of baseline value (0.28 ± 0.09 vs. 0.99 ± 0.14 , $p = 0.01$). Brain–blood flow with 30° HUT was 42% higher compared with the 0° supine position 0.27 ± 0.04 vs. 0.19 ± 0.04 ($p = 0.01$). With 30° HDT brain flow was reduced by 26% to 0.14 ± 0.06 compared to supine 0° values ($p = 0.16$). Heart blood flow was 0.28 ± 0.09 at 0° vs. 0.26 ± 0.06 at 30° HUT ($p = 0.48$) and 0.22 ± 0.07 at 30° HDT ($p = 0.23$) (Fig. 3).

3.3. Blood gas analysis

During L-CPR + ITD the arterial PO_2 was significantly higher with 30° HUT compared with 0° CPR at the same inspired oxygen fraction (175 ± 30 vs. 120 ± 14 , $p = 0.009$) (Table 2). In addition, peak inspiratory pressure was lower with 30° HUT compared with the 0° supine position, 28.7 ± 1.8 mmHg vs. 37.8 ± 1.9 , respectively ($p < 0.001$).




Fig. 3. Blood flow to the heart and brain with supine, 30° head up and 30° head down CPR. CPP: coronary perfusion pressure (CPP), cerebral perfusion pressure (CerPP). * $p < 0.05$ compared with 0°.

3.4. LUCAS alone vs. LUCAS + ITD

To determine the potential importance of the L-CPR + ITD combination vs. L-CPR only, at the end of Protocol A the ITD was removed. We observed an immediate and significant decrease in systolic blood pressure (78 ± 4 vs. 58 ± 7 mmHg, $p = 0.011$), diastolic blood pressure (19 ± 3 vs. 17 ± 2 mmHg, $p = 0.002$), and CPP (25 ± 2 vs. 23 ± 2 mmHg, $p = 0.012$) as soon as the ITD was removed.

3.5. Protocol C

With 0, 10, 20, 30, 40 and 50° HUT, mean ICP decreased linearly as follows: 21 ± 2 , 16 ± 2 , 10 ± 2 , 5 ± 2 , 0 ± 2 , -5 ± 2 respectively, ($p < 0.001$). During the compression phase of CPR the maximum ICP decreased also linearly with 0, 10, 20, 30, 40 and 50° HUT: 30 ± 2 , 24 ± 2 , 16 ± 2 , 12 ± 2 , 7 ± 2 , 4 ± 2 respectively, ($p < 0.001$), whereas CerPP increased linearly ($p = 0.001$), and CPP remained constant (Fig. 4).

Table 2
Respiratory parameters at baseline and between 0, 30° head up and 30° head down CPR.

	Baseline	CPR angle		
		0°	30°	−30°
pH	$7.43 \pm 0.01^*$	7.23 ± 0.02	7.21 ± 0.02	$7.13 \pm 0.03^*$
PCO_2	$41 \pm 1^*$	49 ± 3	46 ± 2	$58 \pm 4^*$
PO_2	104 ± 6	120 ± 14	$175 \pm 30^*$	113 ± 18
HCO_3^-	$27 \pm 1^*$	20 ± 1	$18 \pm 1^*$	$18 \pm 1^*$
BE	$2.9 \pm 0.7^*$	-7.4 ± 0.6	$-9.4 \pm 0.4^*$	$-11.1 \pm 0.6^*$
$ETCO_2$	40 ± 1	35 ± 2	33 ± 2	31 ± 4
ITP max	$3.6 \pm 0.1^*$	7.6 ± 0.7	7.5 ± 0.8	8.7 ± 1.2
ITP min	$2.3 \pm 0.1^*$	-10.3 ± 0.8	-10.3 ± 0.8	$-9.2 \pm 0.9^*$
PIP	$18.7 \pm 0.6^*$	37.8 ± 1.9	$28.7 \pm 1.8^*$	42.7 ± 2.2

ITP: Intrathoracic pressure during chest compression; PIP: peak inspiratory pressure during positive pressure ventilation. Pressures are in mmHg.

* $p < 0.05$ compare to 0° CPR.




Fig. 4. Dose angle relationships with mean aortic pressure and intracranial pressure, coronary perfusion pressure (CPP) and cerebral perfusion pressure (CerPP). *p* value are reported for linear trend.

4. Discussion

An inherent limitation to systemic and cerebral perfusion during conventional CPR is the high resistance generated in the venous circulation with each chest compression, which reduces or eliminates the forward flow pressure gradient to the heart and often to the brain.¹² Results from the current study demonstrate that L-CPR+ITD performed in the head up position (reverse Trendelenberg) significantly increases cerebral perfusion pressure, oxygenation, and cerebral blood flow, compared with either 0° supine or HDT. The shift in the position of the head and body to the upright position immediately reduced ICP and venous pressures, thereby rapidly reducing resistance to forward blood flow generated by L-CPR+ITD. As such, the combination of L-CPR+ITD in the novel HUT position may provide a new way to overcome some of the most fundamental challenges associated with the generation of sufficient circulation to the heart and brain.

During standard CPR in the 0° supine position, arterial and venous pressure rise to a similar degree with each chest compression. This increase in venous pressure is nearly instantaneously transmitted to the brain through venous structures including the paravertebral venous plexus.¹² As shown in Fig. 2, when CPR was performed in the horizontal plane, each compression is associated with a rise in ICP and a concurrent decrease in the aortic minus ICP cerebral perfusion gradient: in some pig studies ICP during the compression phase was >80 mmHg. The findings from the current study, for the first time to our knowledge, demonstrate that tilting the head and body upward during the chest compression and decompression phases reduces venous and ICP pressures thereby significantly increasing brain perfusion and calculated perfusion pressure.

The optimal head up angle is not currently known. From a mechanistic perspective, the results from the current study demonstrate that with each 10° increase in the HUT angle the peak, trough and mean compression phase ICP and right atrial pressure decreased significantly, thus generating a much larger cerebral vascular

pressure gradient between the arteries and veins. The decreases in ICP and right atrial pressure were linear, with a 0.3–0.4 mmHg decrease with each 1° increase in HUT. This resulted in a linear increase in cerebral perfusion pressure. By contrast, HDT resulted in a striking increase in ICP and a decrease in coronary and cerebral perfusion pressures. The absolute ICP values during the compression phase of head down suggest that such an approach may cause significant harm.

In order to consistently perfuse the brain in the elevated head position, it is important to maintain sufficient central blood volume and forward flow. This was accomplished in the current study with L-CPR+ITD, a device combination that harnesses the intrathoracic vacuum generated during the decompression phase of CPR to enhance circulation and further reduce ICP.^{6,10} In this study, without the ITD the systolic and diastolic blood pressure as well as calculated CPP decreased immediately.

This first study on the potential impact of head up position was limited by the experimental design: different interventions were performed in sequential rather than randomized order. In prior studies coronary perfusion pressures have been shown to decrease over time so the absolute increase in these parameters with head-up position may be underestimated.¹³ In addition, this new concept has resulted in a number of currently unanswered questions including: what is the optimal CPR head-up angle, what is the impact of HUT on pulmonary function and circulation, should the entire body be tilted upward or just part of the body, how long can CPR be performed in the head up position, and will survival rates be increased HUT. Studies are currently ongoing to address these important issues.

5. Conclusion

CPR with head up positioning provides a novel way to easily, rapidly, and significantly augment cerebral perfusion during CPR without impairing perfusion to the heart. Mean ICP values

decreased and CerPP increased linearly until 50° head-up positioning. Further studies are warranted to explore this new approach.

Source of funding

This study was funded by an internal grant of the Cardiovascular Division of the University of Minnesota. Guillaume Debaty received a grant from the Region Rhône-Alpes (France) and from the Société Française de Médecine d'Urgence (SFMU) for a post-doctoral fellowship.

Conflicts of interest statement

Keith Lurie is the inventor of the Impedance Threshold Device and is the medical director of Advanced Circulatory System (Roseville, MN). Anja Metzger, and Michael Lick are employees of Advanced Circulatory System.

There are no other conflicts of interest.

Appendix A. Supplementary data

Supplementary data associated with this article can be found, in the online version, at <http://dx.doi.org/10.1016/j.resuscitation.2014.11.019>.

References

- Cohen TJ, Tucker KJ, Lurie KG, et al. Active compression-decompression. A new method of cardiopulmonary resuscitation. Cardiopulmonary Resuscitation Working Group. *JAMA* 1992;267:2916–23.
- Lurie KG, Lindner KH. Recent advances in cardiopulmonary resuscitation. *J Cardiovasc Electrophysiol* 1997;8:584–600.
- Duggal C, Weil MH, Gazmuri RJ, et al. Regional blood flow during closed-chest cardiac resuscitation in rats. *J Appl Physiol* 1993;74:147–52.
- Lick CJ, Aufderheide TP, Niskanen RA, et al. Take heart America: a comprehensive, community-wide, systems-based approach to the treatment of cardiac arrest. *Crit Care Med* 2011;39:26–33.
- Debaty G, Segal N, Matsuura T, et al. Hemodynamic improvement of a LUCAS 2 automated device by addition of an impedance threshold device in a pig model of cardiac arrest. *Resuscitation* 2014. <http://dx.doi.org/10.1016/j.resuscitation.2014.09.013>.
- Metzger AK, Herman M, McKnite S, Tang W, Yannopoulos D. Improved cerebral perfusion pressures and 24-hr neurological survival in a porcine model of cardiac arrest with active compression-decompression cardiopulmonary resuscitation and augmentation of negative intrathoracic pressure. *Crit Care Med* 2012;40:1851–6.
- Lurie KG, Yannopoulos D, McKnite SH, et al. Comparison of a 10-breaths-per-minute versus a 2-breaths-per-minute strategy during cardiopulmonary resuscitation in a porcine model of cardiac arrest. *Respir Care* 2008;53:862–70.
- Ditchey RV, Winkler JV, Rhodes CA. Relative lack of coronary blood flow during closed-chest resuscitation in dogs. *Circulation* 1982;66:297–302.
- Voelckel WG, Lurie KG, Sweeney M, et al. Effects of active compression-decompression cardiopulmonary resuscitation with the inspiratory threshold valve in a young porcine model of cardiac arrest. *Pediatr Res* 2002;51:523–7.
- Lurie KG, Coffeen P, Shultz J, McKnite S, Detloff B, Mulligan K. Improving active compression-decompression cardiopulmonary resuscitation with an inspiratory impedance valve. *Circulation* 1995;91:1629–32.
- Reinhardt CP, Dalhberg S, Tries MA, Marcel R, Leppo JA. Stable labeled microspheres to measure perfusion: validation of a neutron activation assay technique. *Am J Physiol Heart Circ Physiol* 2001;280:H108–16.
- Paradis NA, Martin GB, Goetting MG, et al. Simultaneous aortic, jugular bulb, and right atrial pressures during cardiopulmonary resuscitation in humans. Insights into mechanisms. *Circulation* 1989;80:361–8.
- Lurie KG, Voelckel WG, Zielinski T, et al. Improving standard cardiopulmonary resuscitation with an inspiratory impedance threshold valve in a porcine model of cardiac arrest. *Anesth Analg* 2001;93:649–55.