



EQUIPMENT CARTS

Mobility on castors – complete catalogue



creating medical mobility

CREATING CUSTOMER-SPECIFIC SOLUTIONS:

Creating medical mobility



As one of the world's leading providers of equipment carts and support arms for medical technology, we can look back on more than 20 years of experience in our market, thousands of satisfied customers and countless project implementations.

We focus on the continuous development of our employees' technical and personal skills and on the modular development of our products. This modular system enables us and our customers to build easily their specific products using standard and near-standard elements.





Our slogan "creating medical mobility" aims to make this clear and drive us forward.

Quality you can touch for a lasting cooperation.

Product quality and customer satisfaction are of paramount importance to ITD GmbH. In order to fulfill these high quality expectations, ITD GmbH is committed to the most modern technologies.

The combination of knowing the fields of application for the products, market-oriented development and above all the close partnership with the customer - from idea to finished product - a quality product is created to meet individual requirements. Today as well as in future.

CONTENTS

	Product attributes and applications	4
	Overview of product portfolio	6
	uni-cart	10
	Application examples uni-cart	12
	Order configurations uni-cart	18
	System components uni-cart.....	20
	vexio-cart	30
	Application examples vexio-cart	32
	Order configurations vexio-cart	38
	System components vexio-cart	40
	pro-cart	58
	Application examples pro-cart	60
	Order configurations pro-cart ..	66
	System components pro-cart	68
	classic-cart	78
	Application examples classic-cart	80
	Order configurations classic-cart	86
	System components classic-cart	90
	compact-cart	102
	Application examples compact-cart	104
	Order configurations compact-cart	110
	System components compact-cart	112



Monitor carts 124

Order configurations	126
----------------------------	-----

Customised product solutions cart-couture 133

Accessories 136



Pivot arms rm-port	140
Variable height support arms flexion-port	141-143
Monitor crossbeam / Monitor shelf	144
Counterweight set	145
Baskets	146
Computer holder	147
Laptop holder	148
Power pack holder	148
Adapter plates for universal adapter	149
Endoscope holder	150
Drip glass	151
Connector holder	151
Rinsing bottle holder	152
Hose holder	152
Camera holder / Universal holder	152
(Gel)-bottle holder	153
Gas bottle holder	154
IV pole	155
Standard runner	155
Anchor / clamps	156
Rail fastener	156
Holder for suction system	157
Adapter pieces	157
Foot switch cradle	158
Jointed arms	158
Cable winders	159
Cable holder	159
Connection cable	160
EURO cable	161
Adapter cable	161
Equipotential plug POAG	162
Equipotential cable POAG	162
Socket strip attachment	163
Plug removal barrier	164
Spring nut	165
Equipment fixing corners	165

Item Number Index 166

Keyword Index 170

OUR EXPERIENCE – YOUR SUCCESS

The objective of ITD is to react promptly to the latest market and environmental requirements. The aim is to allow our customers to exploit new opportunities for success. This is why we are dedicated to channelling our experience and our “craft assets” – innovative excellence, technical expertise and design competence into our products. ITD creates new benchmarks with our products feature attributes like modularity, ergonomic design and quality. They are supported by service and customizing. This powerful combination generates product solutions that meet all the product requirements of our sophisticated customers.


Modularity – smart and individual

- The standard building-block system permits individual and cost-effective configurations.
- The modular support systems can be upgraded at any time. This permits prompt reactions to new demands in the marketplace and guarantees future-proof investments.
- The variability of height and / or width combine with an extensive range of generic accessories overarching all product lines to meet diverse needs and requirements.
- Easy planning is due to online product configurator everytime possible.



Ergonomic design – adaptable and natural

- Smooth-running twin castors ensure simple manoeuvring of the equipment carts as well as the safe rolling over of obstacles.
- Adapted working areas and handle positions allow users to maintain a natural body position (sit / stand workstation).
- The support systems have been adapted to natural sequences of movement. The design of the base provides adequate legroom and permits freedom of movement.



Hygiene and Design – focused and modern

- The link between design and functionality focuses on instrumentation.
- The modern look and appealing, streamlined design underscore the high quality and value of the instruments.
- The absolute high standards of hygiene required in the everyday routine of hospitals and medical practices are reflected in the detail of the practical design.





Quality – long service life and safety

- The use of high-quality materials, such as zinc-plated and powder-coated sheet steel or aluminium extrusions, guarantees the long service life and dependability of the product solutions.
- All support systems comply with the European Medical Device Regulation (MDR) 2017/745 as well as the standards according to DIN EN 60601-1 (3rd Edition).
- Quality assurance is carried out strictly in conformity with DIN EN ISO 13485. Electrical safety is also tested and documented down to the last detail.



Service – personal and efficient

- Professional technical personnel provide friendly and expert advice. Personal consultations are given on site or at our showrooms.
- Support is delivered throughout the project period by a dedicated personal contact.
- Expert foreign agents provide the foundation for the international presence of ITD.
- Detailed offers are prepared promptly.
- National and international logistics are supported comprehensively by e. g. ex-stock items and deliveries to end customers.



Customizing – specific and characteristic

- Tailor-made customised solutions are based on an extensive standard portfolio enabling fast response to requirements at economical prices.
- Design of the equipment can be made in the colour of choice or with the corporate identity as required by the customer.
- On request, a package can be provided with customised corporate and product labelling, e.g. with the company logo or barcode.
- The design of the individual equipment carriers is carried out in cooperation with the ITD representative.

OUR EQUIPMENT CARTS



Surgery, intensive care, anaesthesia or laboratory settings – each application presents different challenges to an equipment cart. Excellent that ITD can offer an optimum solution for each scenario. The ITD product portfolio is convincing point by point through versatility and sophistication.

The following tabular presentation provides you with a quick initial overview of possible equipment carts and can assist you in making the right preliminary decision.



classic-cart



uni-cart

AN OVERVIEW OF OUR EQUIPMENT CARTS



■ uni-cart



■ vexio-cart

Load capacity	with 4 stabilisers	with 5 stabilisers	
Total load capacity:	65 kg / 143 lbs	50 kg / 110 lbs	65 kg / 143 lbs
Load capacity of shelf:	10 kg / 22 lbs distributed load		15 kg / 33 lbs distributed load
Dimensions - Basic Frame			
Width in mm:	540	670	530
Depth in mm:	580	650	610
Total height in mm:	820 / 1145 / 1305		830 / 998 / 1160 / 1320
Dimensions - Shelf			
Width in mm:	400		424
Depth in mm:	360		430
Electrification			
Isolating transformer:	660 VA		660 VA
Central on / off switch:	no		yes
Socket strips:	yes		yes
Ordering data			
Catalogue page:	18 - 29		38 - 56



■ pro-cart



■ classic-cart



■ compact-cart

80 kg / 176 lbs

20 kg / 44 lbs distributed load

150 kg / 330 lbs

50 kg / 110 lbs distributed load

180 kg / 396 lbs

50 kg / 110 lbs distributed load

610

620

884 / 1172 / 1332

614 / 714 / 1214

625

829 / 1117 / 1437

618 / 668 / 818

624

1146 / 1306 / 1466

450

400

450 / 500 / 600

491

450 / 500 / 650

491

660 VA / 1200 VA

yes

yes

660 VA / 1000 VA / 1600 VA

yes

yes

1200 VA / 1600 VA / 2000 VA

yes

yes

66 - 77

86 - 101

110 - 123



Support arms - Flexibility in the room complete catalogue

You will find information on our support arms in the dedicated separate catalogue. You can request your copy of this catalogue from the address below or download it from the internet on www.itd-cart.com.

ITD GmbH
Sportplatzstraße 3
84381 Johanniskirchen
Germany

Fon: +49 89 61 44 25-0
Fax: +49 89 61 44 25-200

www.itd-cart.com
sales@itd-cart.com