

The EB-J10 Scopes

Experience crystal clarity – High-Definition at its best



The PENTAX Medical Video Bronchoscope EB-J10 product range offers sharp and clear High-Definition image quality. Combined with a High-Definition processor offering i-scan technology, the bronchoscopes support faster detection, precise demarcation and support the adequate patient treatment.

The unmatched ergonomic design and flexibility of the bronchoscopes offer you ease of operation in daily use. Complying with increasing hygiene requirements and STERRAD® compatibility the EB-J10 bronchoscopes provide patient and user safety.

A ground-breaking
endoscopic
advancement



Sharp imaging for detailed visualization

The excellent visualization of HD offers fast and detailed bronchial system investigation and leads to more precise detection of abnormalities supporting precise demarcation of lesions or tumors.

Reliable for user and patient safety

For user and patient safety, the EB-J10 bronchoscopes comply with increasing hygiene requirements through single-use accessories and STERRAD® compatibility.

Unmatched ergonomic design for ease of operation

Unmatched ergonomics and lightness of the control body in combination with the excellent scope maneuverability thanks to a very small distal tip bending radius ease operation.

Video Bronchoscopes

	EB15-J10	EB19-J10
Distal end diameter [mm]	5.4	6.1
Insertion tube diameter [mm]	5.2	6.4
Instrument channel diameter [mm]	2.0	2.8
Working length [mm]	600	600
Field of view [°]	120	120
Depth of field [mm]	3-100	3-100
Angulation up/down [°]	210/130	180/130

EC REP EMEA Headquarter
Germany

PENTAX Europe GmbH
Julius-Vosseler-Straße 104 . 22527 Hamburg
Tel.: +49 40 / 5 61 92 - 0 . Fax: +49 40 / 5 60 42 13
E-mail: info.emea@pentaxmedical.com . www.pentaxmedical.com

TÜV Süd CE0123 • Medical device class: IIa • This product must be used only by Healthcare professionals. Before usage and for detailed product specifications, please refer to the instructions for use. In the interest of technical process, specifications may change without notice.

PENTAX
MEDICAL
Excellence in Focus