



EASY PULSE®

Efficient and compact
mechanical CPR thanks to
technological innovation



SCHILLER
The Art of Diagnostics

EASY PULSE®

SCHILLER has taken the next step in improving CPR quality further with the launch of the smallest mechanical chest compression device: the EASY PULSE.

Performing manual chest compressions well for an extended period of time is almost impossible. Not only is it physically demanding, but other actions are also required simultaneously, such as vital signs monitoring.

SCHILLER's EASY PULSE is the solution for more effective resuscitation: this portable, stand-alone device delivers chest compressions automatically at a consistent rate and depth. It is directly attached to the patient's upper body and can thus be used in any situation, regardless of ambient conditions.

The EASY PULSE is such a simple device that a single person can operate it – whether a paramedic or a lay rescuer. Its size and weight are unparalleled: the EASY PULSE weighs less than 3.5 kg, opening up entirely new possibilities e.g. for air rescue services.



SLIDER

The slider makes it possible for a single person to prepare the patient for mechanical chest compressions. It also takes very little time to attach the device thanks to the practical buckle system - less than 15 seconds for trained personnel.



CHARGER

The battery can be charged in a dual charger (less than 2 hours) or in turbo mode (less than 50 minutes). The EASY PULSE is powered by a rechargeable battery (running time approx. 45 minutes).



Technological innovations

Lightweight

- World's smallest and lightest electrically driven CPR device, weighing less than 3.5 kg
- Ideal for limited spaces

Comprehensive

- Multidirectional chest compressions
- High-quality chest compressions of consistent rate and depth, which is impossible to achieve with manual compressions
- Possibility of 30:2 compression - ventilation cycle

Simple

- Thanks to slider and buckle system, a single person can easily attach the device
- Easy to position - easy to operate

Additional features

- Replaceable battery with charge level indicator
- Autonomous operation for 45 minutes (with one battery)
- Connections: external DC input, USB



EXTERNAL CONNECTION CABLE

The EASY PULSE can be directly operated using a 12VDC - 24VDC power source via a special connection cable. A power supply unit for 100VAC to 240VAC input voltage is supplied.

EASE OF USE

In an emergency every move has to be right. The EASY PULSE is simple and intuitive to operate: pressing a button twice is enough to initiate resuscitation immediately.





Suomi – Finland
Medifin Healthcare Oy
Puh. +358 (0)10 321 3870
asiakaspalvelu@medifin.fi
www.medifin.fi



Asia
SCHILLER Asia-Pacific / Malaysia
52200 Kuala Lumpur, Malaysia
Phone +603 6272 3033
Fax +603 6272 2030
sales@schiller.com.my
www.schiller-asia.com



Austria
SCHILLER Handelsgesellschaft m.b.H.
A-4040 Linz
Phone +43 732 709 90
Fax +43 732 757 000
sales@schiller.at
www.schiller.at



China
Alfred Schiller Medical Equipment Co. Ltd.
100015 Beijing, China
Phone +86-10-84565453
info@schillermedical.com
www.schillermedical.cn



France
SCHILLER Médical S.A.S.
F-67162 Wissembourg/Cedex
Phone +33 3 88 63 30 00
Fax +33 3 88 94 12 82
info@schiller.fr
www.schiller-medical.com



France (distribution France)
SCHILLER France S.A.S.
F-77600 Bussy St Georges
Phone +33 1 64 66 50 00
Fax +33 1 64 66 50 10
infoschiller@schiller-france.fr
www.schiller-france.com



Germany
SCHILLER Medizintechnik GmbH
D-85622 Feldkirchen b. München
Phone +49 89 62 99 81-0
Fax +49 89 62 99 81-54
info@schillermed.de
www.schillermed.de



Hungary
SCHILLER Diamed Ltd.
H-1141 Budapest
Phone +36 (1) 383-4780 / 460-9491
Fax +36 (1) 383-4778
sales@schiller.at
www.schiller-hungary.hu



India
SCHILLER Healthcare India Pvt. Ltd.
Mumbai - 400 001, India
Phone +91 22 6152 3333/ 2920 9141
Fax +91 22 2920 9142
sales@schillerindia.com
www.schillerindia.com



Croatia
Schiller medicinski instrumenti d.o.o.
10000 Zagreb
Phone +385 1 309 66 59
Fax +385 1 309 66 60
info@schillerz.hr
www.schiller.ch



Latin America
SCHILLER Latin America, Inc.
Doral, Florida 33172
Phone +1 954 673 0358
Fax +1 786 845 06 02
info@schillerla.ch
www.schillerla.ch



Poland
SCHILLER Poland Sp. z o.o.
PL-02-729 Warszawa
Phone +48 22 8432089
Fax +48 22 8432089
schiller@schiller.pl
www.schiller.pl



Russia & C.I.S.
SCHILLER AG Rep. office
125124 Moscow, Russia
Phone +7 (495) 970 11 33
Fax +7 (495) 970 11 33
mail@schiller-ag.com
www.schiller-cis.com



Serbia
SCHILLER d.o.o.
11010 Beograd
Phone +381 11 39 79 508
Fax +381 11 39 79 518
info@schiller.rs
www.schiller.rs



Slovenia
SCHILLER d.o.o.
2310 Slovenska Bistrica
Phone +386 2 843 00 56
Fax +386 2 843 00 57
info@schiller.si
www.schiller.si



Spain
SCHILLER ESPAÑA, S.A.
E-28230-Las Rozas/Madrid
Phone +34 91 713 01 76
Fax +34 91 355 79 33
schiller@schiller.es
www.schiller.es



Switzerland
SCHILLER-Reomed AG
CH-8953 Dietikon
Phone +41 44 744 30 00
Fax +41 44 740 37 10
sales@schiller-reomed.ch
www.schiller-reomed.ch



Turkey
SCHILLER TÜRKİYE
Okmeydanı-Sisli - İstanbul
Phone +90 212 210 8681 (pbx)
Fax +90 212 210 8684
sales@schiller-turkiye.com
www.schiller-turkiye.com



USA
SCHILLER America Inc.
Doral, Florida 33172
Phone +1 786 845 0620
Fax +1 786 845 06 02
sales@schilleramerica.com
www.schilleramerica.com



SCHILLER
The Art of Diagnostics