


ERG 910 PLUS ERG 911 PLUS

Functional, high-quality ergometers



SCHILLER
The Art of Diagnostics

ERG 910 PLUS

ERG 911 PLUS

- ∴ Pivoting display 180°
- ∴ Easy to operate
- ∴ Almost noiseless operation
- ∴ Low access
- ∴ Stable steel tube construction



APPLICATION AREAS

The bike ergometers ERG 910 plus and ERG 911 plus were developed for use in the cardiology sector based on the proven previous models, with particular care to ensure highest stability and practicability. The devices are used for performing defined exercises during ergometric examinations and therapeutic applications.

ERG 911 plus is available in both versions "easy" and "smart". "Easy" devices are intended to be used with a master device (remote operation: ECG, PC, etc.) whereas "smart" models can be used as stand-alone devices. Furthermore, many options can be added to the available features on both devices.

Regardless of their features, the devices meet the highest quality standards in compliance with DIN VDE 0750-238 for accurate physical exertion tests to conduct measurements in cardiovascular and pulmonary function diagnostics.

- ∴ Blood pressure measurement using QRS trigger
- ∴ Electrical saddle height adjustment
- ∴ SpO₂ measurement
- ∴ Bluetooth connection
- ∴ ECG amplifier module



ERG 910 PLUS – EASY

- ❖ Remote operation (type "easy")
- ❖ load range 20-500 Watts
- ❖ optional: power supply 110-120V
- ❖ optional: 20-999 Watt



ERG 911 PLUS – EASY / SMART

- ❖ remote (type "easy") and stand-alone (type "smart") operation
- ❖ load range 20-999 Watts
- ❖ various optional extensions available for both models (see technical data)



EASY

Ergometers equipped with an "easy" control panel, on which they receive exercise specifications from a master device (ECG, PC, ...), are used for remote operation. All major exercise data (such as blood pressure, heart rate, pedalling rate, instruction concerning correct pedalling rate) is displayed in a compact format. The optional features, such as blood pressure measurement and electric saddle height adjustment, can also be conveniently adjusted from the Touch display.



SMART

Ergometers equipped with a "smart" control panel, can be used for remote as well as for stand-alone operation. Exercise instructions can be transmitted from a master device (ECG, PC, ...) or from a programmable ergometry program located on the ergometer. Further, training programs with exercise types such as pulse steady state, constant load and interval training can also be programmed.



Americas
SCHILLER Americas Inc.
 Doral, Florida 33172

North America:
 Phone +1 786 845 06 20
 Fax +1 786 845 06 02
sales@schilleramerica.com
www.schilleramerica.com

Latin America & Caribbean
 Phone +1 305 591 11 21
 Fax +1 786 845 06 02
sales@schillerla.ch
www.schillerla.ch



Asia
SCHILLER Asia-Pacific / Malaysia
 52200 Kuala Lumpur, Malaysia
 Phone +603 6272 3033
 Fax +603 6272 2030
sales@schiller.com.my
www.schiller-asia.com



Austria
SCHILLER Handelsgesellschaft m.b.H.
 A-4040 Linz
 Phone +43 732 70 99
 Fax +43 732 757 000
sales@schiller.at
www.schiller.at



China
Alfred Schiller (Beijing) Medical Equipment Co.,Ltd.
 100015 Beijing, China
 Phone +86 010 52007020
 Fax +86 010 52007020-8016
info@schillerchina.com
www.schillermedical.cn



Croatia
SCHILLER d.o.o.
 10000 Zagreb
 Phone +385 1 309 66 59
 Fax +385 1 309 66 60
info@schillerzg.hr
www.schiller.hr



France
SCHILLER Médical
 F-67160 Wissembourg
 Phone +33 3 88 63 36 00
 Fax +33 3 88 63 36 49
info@schiller.fr
www.schiller-medical.com



France (distribution France)
SCHILLER France S.A.S.
 F-77608 Bussy St Georges
 Phone +33 1 64 66 50 00
 Fax +33 1 64 66 50 10
contact@schillerfrance.fr
www.schiller-france.com



Germany
SCHILLER Medizintechnik GmbH
 D-85622 Feldkirchen b. München
 Phone +49 89 62 99 81 0
 Fax +49 89 609 50 90
info@schillermed.de
www.schillermed.de



Hungary
SCHILLER Diamed Ltd.
 H-1141 Budapest
 Phone +36 1 383 4780 / +36 1 460 9491
 Fax +36 1 383 4778
info@schillerhungary.hu
www.schillerhungary.hu



India
SCHILLER Healthcare India Pvt. Ltd.
 Mumbai – 400 059, India
 Phone +91 22 6152 3333 / 2920 9141
 Fax +91 22 2920 9142
sales@schillerindia.com
www.schillerindia.com



Poland
SCHILLER Poland Sp. z o.o.
 PL-02-729 Warszawa
 Phone +48 22 843 20 89 / +48 22 647 35 90
 Fax +48 22 843 20 89
schiller@schiller.pl
www.schiller.pl



Russia & C.I.S.
AO SCHILLER.RU
 125040 Moscow, Russia
 Phone +7 (495) 970 11 33
 Fax +7 (495) 956 29 10
mail@schiller.ru
www.schiller.ru



Suomi – Finland
Medifin Healthcare Oy
Puh. +358 (0)10 321 3870
asiakaspalvelu@medifin.fi
www.medifin.fi



Serbia
SCHILLER d.o.o.
 11010 Beograd
 Phone +381 11 39 79 508
 Fax +381 11 39 79 518
info@schiller.rs
www.schiller.rs



Slovenia
SCHILLER d.o.o.
 2310 Slovenska Bistrica
 Phone +386 2 843 00 56
 Fax +386 2 843 00 57
info@schiller.si
www.schiller.si



Spain
SCHILLER ESPAÑA, S.A.
 ES-28232 Las Rozas/Madrid
 Phone +34 91 713 01 76
 Fax +34 91 355 79 33
schiller@schiller.es
www.schiller.es



Switzerland
SCHILLER-Reomed AG
 CH-8953 Dietikon
 Phone +41 44 744 30 00
 Fax +41 44 740 37 10
sales@schiller-reomed.ch
www.schiller-reomed.ch



Turkey
SCHILLER TÜRKİYE
 Okmeydanı-Sisli – Istanbul
 Phone +90 212 210 8681 (pbx)
 Fax +90 212 210 8684
info@schiller.com.tr
www.schiller-turkiye.com



United Kingdom
SCHILLER UK
 Bellshill, ML4 3PR
 Phone +44 1698 744 505
 Fax +44 1698 744 474
sales@schilleruk.com
www.schilleruk.com

Article no.: 2.500772
 2018/03_EN



SCHILLER
 The Art of Diagnostics