



CARDIOVIT CS-200 Office

Wireless ECG acquisition




CARDIOVIT CS-200 OFFICE

The CARDIOVIT CS-200 Office offers cutting-edge diagnostic capability for 12-lead resting and stress wireless ECG acquisition. The choice is yours: Whether only the software for any laptop or as a package integrated in a medical PC, it's the perfect solution for GPs, clinics and hospitals.



MS-12 BLUE



**SWITCH BETWEEN TWO CUSTOMISABLE
SCREEN LAYOUTS AT THE TOUCH OF A
BUTTON**

ECG ACQUISITION WITH MAXIMUM FREEDOM

With the introduction of the MS-12 blue, SCHILLER provides a new solution that greatly simplifies ECG acquisition.

- ❖ Allows to freely operate in any environment without worrying about cable length.
- ❖ Weighing only 115g and fitting in one hand, the MS-12 blue sits on the bed or on the patient body during ECG acquisition.
- ❖ Motion artefact is minimized without an analogue cable connection and an extremely lightweight patient cable.
- ❖ One battery charge guarantees continuous operation for more than 36 hours and more than a week with the power-saving mode.

ADVANCED CLINICAL TOOLS

The CS-200 Office can be supplemented with high-end software features:

- ❖ ETM interpretation
- ❖ ETM sport
- ❖ 12-lead exercise ECG
- ❖ 3D vector ECG
- ❖ Arrhythmia detection and classification
- ❖ Risk factors

USER-DEFINED WORKFLOW (WITH PRINTOUTS OR PAPERLESS)

The CARDIOVIT CS-200 Office fulfills the needs of every user thanks to a customizable user interface, stage and report printout with laser printer, and seamless connectivity to EMR, PACS, HIS or SCHILLER SEMA data management system.

SOLUTIONS FOR ANY ENVIRONMENT

The CS-200 Office can be easily installed on your laptop or PC, as well as purchased as a ready-to-use system including a robust trolley and a 21.5" touch panel PC with pre-installed software (see picture).

Note:

CARDIOVIT CS-200 Office is also available with MS-12 USB (USB connection).

High-resolution and easy-to-clean
21.5" touch screen.

Intuitive and clear user interface
suitable for touch screen.

Soft handle positioned at an
adjustable height.

Support for PC power
supply and cables.

Transparent bin to store and
transport medical items.

User-defined screen layouts
and workflow.

MS-12 blue: Small but
powerful wireless ECG
acquisition module with
on-screen signal view and
hookup advisor.



Easy-glide wheels to smoothly
maneuver the trolley.



Americas
SCHILLER Americas Inc.
 Doral, Florida 33172

North America:
 Phone +1 786 845 06 20
 Fax +1 786 845 06 02
sales@schilleramerica.com
www.schilleramerica.com

Latin America & Caribbean
 Phone +1 305 591 11 21
 Fax +1 786 845 06 02
sales@schillerla.ch
www.schillerla.ch



Asia
SCHILLER Asia-Pacific / Malaysia
 52200 Kuala Lumpur, Malaysia
 Phone +603 6272 3033
 Fax +603 6272 2030
sales@schiller.com.my
www.schiller-asia.com



Austria
SCHILLER Handelsgesellschaft m.b.H.
 A-4040 Linz
 Phone +43 732 70 99
 Fax +43 732 757 000
sales@schiller.at
www.schiller.at



China
Alfred Schiller (Beijing) Medical Equipment Co., Ltd.
 100015 Beijing, China
 Phone +86 010 52007020
 Fax +86 010 52007020-8016
info@schillerchina.com
www.schillermedical.cn



Croatia
SCHILLER d.o.o.
 10000 Zagreb
 Phone +385 1 309 66 59
 Fax +385 1 309 66 60
info@schillerzg.hr
www.schiller.hr



France
SCHILLER Médical S.A.S.
 F-67162 Wissembourg
 Phone +33 3 88 63 36 00
 Fax +33 3 88 63 36 49
info@schiller.fr
www.schiller-medical.com



France (distribution France)
SCHILLER France S.A.S.
 F-77608 Bussy St Georges
 Phone +33 1 64 66 50 00
 Fax +33 1 64 66 50 10
contact@schillerfrance.fr
www.schiller-france.com



Germany
SCHILLER Medizintechnik GmbH
 D-85622 Feldkirchen b. München
 Phone +49 89 62 99 81 0
 Fax +49 89 609 50 90
info@schillermed.de
www.schillermed.de



Hungary
SCHILLER Diamed Ltd.
 H-1141 Budapest
 Phone +36 1 383 4780 / +36 1 460 9491
 Fax +36 1 383 4778
info@schillerhungary.hu
www.schillerhungary.hu



India
SCHILLER Healthcare India Pvt. Ltd.
 Mumbai - 400 059, India
 Phone +91 22 6152 3333 / 2920 9141
 Fax +91 22 2920 9142
sales@schillerindia.com
www.schillerindia.com



Poland
SCHILLER Poland Sp. z o.o.
 PL-02-729 Warszawa
 Phone +48 22 843 20 89 / +48 22 647 35 90
 Fax +48 22 843 20 89
schiller@schiller.pl
www.schiller.pl



Russia & C.I.S.
AO SCHILLER.RU
 125124 Moscow, Russia
 Phone +7 (495) 970 11 33
 Fax +7 (495) 956 29 10
mail@schiller.ru
www.schiller.ru
www.schiller-cis.com



Suomi - Finland
Medifin Healthcare Oy
 Puh. +358 (0)10 321 3870
asiakaspalvelu@medifin.fi
www.medifin.fi



Serbia
SCHILLER d.o.o.
 11010 Beograd
 Phone +381 11 39 79 508
 Fax +381 11 39 79 518
info@schiller.rs
www.schiller.rs



Slovenia
SCHILLER d.o.o.
 2310 Slovenska Bistrica
 Phone +386 2 843 00 56
 Fax +386 2 843 00 57
info@schiller.si
www.schiller.si



Spain
SCHILLER ESPAÑA, S.A.
 ES-28232 Las Rozas/Madrid
 Phone +34 91 713 01 76
 Fax +34 91 355 79 33
schiller@schiller.es
www.schiller.es



Switzerland
SCHILLER-Reomed AG
 CH-8953 Dietikon
 Phone +41 44 744 30 00
 Fax +41 44 740 37 10
sales@schiller-reomed.ch
www.schiller-reomed.ch



Turkey
SCHILLER TÜRKİYE
 Okmeydanı-Sisli - Istanbul
 Phone +90 212 210 8681 (pbx)
 Fax +90 212 210 8684
info@schiller.com.tr
www.schiller-turkiye.com



United Kingdom
SCHILLER UK
 Bellshill, ML4 3PR
 Phone +44 1698 744 505
 Fax +44 1698 744 474
sales@schilleruk.com
www.schilleruk.com

